



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Product Overview

The Eagle 3000 Series Sterilizers feature an advanced microcomputer control system, providing the latest standard for cycle setup, selection and monitoring. Once settings are made and the cycle started, microcomputers accurately monitor and control system operations. The Control Column includes an easy-to-read vacuum fluorescent display and an integral thermal printer.

Eagle 3000 Series Sterilizers can be used in a hospital or laboratory and offer a choice of two configurations:

Gravity

For sterilizing liquids, in flasks with vented closures, at 212° to 254°F (100° to 122°C), and also heat and moisture stable goods at 212° to 280°F (100° to 138°C).

Prevacuum

For efficient, high-volume processing of heat-and moisture-stable materials, such as fabrics, wrapped hard goods, glassware or animal bedding at 270° to 280°F (132° - 138°C) through prevacuum cycle. Prevacuum also provides gravity-type cycles.



**AMSCO EAGLE 3000 SERIES
MEDIUM STERILIZERS**

Specifications:

Single Door

Product Configurations:

3031	24 x 36 x 36	Gravity
3041	24 x 36 x 48	Gravity
3051	24 x 36 x 60	Gravity
3033	24 x 24 x 36	Prevacuum
3043	24 x 36 x 48	Prevacuum
3053	24 x 36 x 60	Prevacuum

Double Door

Product Configurations:

3041	24 x 36 x 48	Gravity
3051	24 x 36 x 60	Gravity
3043	24 x 36 x 48	Prevacuum
3053	24 x 36 x 60	Prevacuum

Door Hinge

On Right Side
On Left Side

Mounting

Cabinet Enclosed/Freestanding
Recessed

Double Door Mounting

Recessed through One Wall
Recessed through Two Walls