



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Product Overview

The most accurate delivery of anesthetic agents over a flow range of 300 mls to 15L, which provides excellent reliability.

Options:

- Forane 19.1 vaporizer
- Halothane 19.1 vaporizer
- Ethrane 19.1 vaporizer



DRAGER 19.1 VAPORIZERS

Specifications:

Weight:.....6 kg
Height:.....205 mm
Diameter:.....100 mm
Temperature range:.....10 - 40 degrees Celsius
Flow Range:.....0.3 - 15l/min

Pressure Compensation Range:

Pressure difference at pressure change:.....20 mbar
Maximum Pressure:.....200 mbar
Maximum Inclination:.....45 degrees