



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Product Overview

The Narkomed 2B offers unsurpassed quality and reliability in gas delivery and ventilation, coupled with built-in capabilities for external monitoring. Narkomed 2B is a dependable, user-friendly foundation upon which to build a multi-faceted anesthesia workstation.

Incorporating ergonomic design principles, Narkomed 2B appeals to a broad range of users. Monitor functions and alarm limits provide optimal control, and are complemented by Drager's central alarm network. Flexible design and a wide range of options enable the Narkomed 2B to meet custom user needs in an extensive variety of applications, and provide hospitals with long term reliability and reduced service costs.



DRAGER NARKOMED 2B

Specifications:

Standard Features Include:

- AV-E Electronic Ventilator
- Dual Flowmeters
- Flow Control Knob Guard
- Touch Coded Oxygen Flow Control Knob
- DISS Pipeline Inlets
- Pipeline Pressure Gauges
- Cylinder Yoke for Oxygen
- Cylinder Yoke for Nitrous Oxide
- Cylinder Pressure Gauges
- Triple Vaporizer Exclusion Valve
- Fresh Gas Outlet with Locking Device
- Minimum Oxygen Flow System
- Oxygen Flush Control
- Oxygen Failure Protection Device
- Oxygen Ratio Controller
- Flowmeter Lights
- 3 Drawers with Writing Tray and Central Locking System
- Footrest, Clipboard Hook and Cylinder Cradles