



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

**Product Overview**

A modern anesthetic designed for “out-of-circuit” use in continuous flow techniques of inhalation anesthesia.

Equipped in temperature-compensated and pressure-compensated capabilities.

**Options:**

- Forane
- Halothane
- Ethrane



**OHMEDA TEC 4 VAPORIZER**

**Specifications:**

**Flow range (LPM):**.....0.25-15  
**Ambient temperature range:**.....18-35 degrees Celsius  
**Capacity (dry wick):**.....135 cc  
**Capacity (wet wick):**.....100

**Calibration ranges:**

Enflurane:.....0-7%  
 Halothane:.....0-5%  
 Isoflurane:..... 0-5%

**Adjustments increments:**     .25 (0-0.5) and .5 (0.5-5)  
**Allowable tilt:**.....180 degrees