



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Product Overview

It offers significant performance and convenience benefits, coupled with important safety features, to improve the clinical management of volatile anesthetic agents in continuous flow techniques.

Options:
Forane
Ethrane
Halothane
Sevoflurane



OHMEDA TEC 5 VAPORIZER SPECIFICATIONS

Calibration: The dials of all Tec 5 vaporizers are calibrated to 5% v/v in steps of 0.2% up to 1%, and in steps of 0.5% from 1% to 5%

Resistance to gas flow
"Off" Isolated - no resistance from vaporizer
"On" 5 L/min O₂

Operating temperature
18 - 35 degrees Celsius

Liquid capacity
To charge dry vaporizers: 300 mL (nominal)
Retained by wick system: 75 mL (nominal)

Weight and Dimensions:
Height: 237 mm
Width: 114 mm
Depth: 97 mm
Weight: 7.0 kg