



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Product Overview

Designed specifically for the delivery of Desflurane anesthetic agent.

Heats Desflurane to maintain temperature and vapor pressure for consistent output.



OHMEDA TEC 6 VAPORIZER

Specifications:

Height:.....25 cm (9.9")
Depth:.....23.5 cm (9.3")
Width:.....11 cm (4.3")
Weight:.....9.5 kg (20 lbs)

Concentration range:.....1 to 18%
Calibrated flow range:.....0.2 to 10 L/min (at 21 degrees Celsius)
Max. permissible flow:..... 20 L/min

Liquid capacity:.....375 mL working volume
.....250 mL refill indicator
..... 50 mL reserve

Power consumption:..... 250 V A nominal
Voltage requirement:...either 220 to 240 V or 100 to 120 V a.c 50/60 Hz
Battery type:.....Duracel MN/1604