



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Microscope Overview

Notable Features of the OPMI-6 are:

- Motorized Focus and Zoom
- New Fiber-Optic Light Source and Cable
- Foot pedal for Motorized Features
- Various Optics Available – OPMI 6, 6CFC, 6S

Spare Bulbs, Sterilizable Knob Covers and Dust Cover provided with Every Microscope System.

2 Year Warranty Including Parts and Labor.



ZEISS OPMI - 6 SERIES

Specifications:

Binocular:.....Inclined Inclinable Upgrade Available

Objective Lens:.....Your Choice: 150mm-400mm

Magnification :..... 1:4 Ratio; Motorized

Focused:.....Motorized, 50mm Range

Eyepieces:.....12.5X HEP

Pupillary Distance:.....50mm-70mm

Filter:.....Standard UV

Light Source:.....Variable 15V, 150W

Illuminated Field:.....Up to 48mm

Power Source:..... 110V-240V; 50-60 Hz

Mounting System:..... Floorstand w/Locking Castors

Horizontal Reach:..... Arm 24 inches

End of Optics:..... 29 inches

Vertical Range (Macro Focus):.....65 inches

Arm:.....Pole Style, Sliding

Shipping Weight:.....200 lbs

Warranty Term:.....2 Years Parts and Labor