



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Microscope Overview

Notable Features of the OPMI®-1 are:

- Exceptional Clarity and Resolution
- Five Steps of Magnification and Fine Focus
- Rollable Floor Stand
- Cool, Fiber Optic Illumination

Spare Bulbs, Sterilizable Knob Covers and Dust Cover provided with Every Microscope System.

2 Year Warranty Including Parts and Labor.



ZEISS OPMI® - 1

Specifications:

Microscope Body:.....OPMI® 1-FC

Magnification:.....2.7, 4.1x, 6.8x, 10.9x, 17x

Fine Focus:.....40mm Excursion

Binoculars & Eyepieces:

Binoculars are straight

F=170

Eyepieces are 10x/22B widefield,
High eyepoint with adjustable eyecups.

Objective:.....f=250mm

Couplings:.....Inclined coupling

Illumination:.....12V 100W fiber optic illumination

Suspension:.....S-1 Floorstand

Weight Restriction:.....7 kg maximum

System Height:.....1677 mm (65.91 inches)
For Floorstand

Warranty Term:.....2 Years Parts and Labor

Optional Accessories:

Co-Observation:.....Beam splitter with monocular
observation available