



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

### Microscope Overview

Notable Features of the OPMI®-1 are:

- Exceptional Clarity and Resolution
- Five Steps of Magnification and Fine Focus
- Rollable Floor Stand
- Cool, Fiber Optic Illumination

Spare Bulbs, Sterilizable Knob Covers and Dust Cover provided with Every Microscope System.

2 Year Warranty Including Parts and Labor.



**ZEISS OPMI® - 1**

### Specifications:

**Microscope Body:**.....OPMI® 1-FC

**Magnification:**.....2.7, 4.1x, 6.8x, 10.9x, 17x

**Fine Focus:**.....40mm Excursion

#### Binoculars & Eyepieces:

Binoculars are straight

F=170

Eyepieces are 10x/22B widefield,  
High eyepoint with adjustable eyecups.

**Objective:**.....f=250mm

**Couplings:**.....Inclined coupling

**Illumination:**.....12V 100W fiber optic illumination

**Suspension:**.....S-1 Floorstand

**Weight Restriction:**.....7 kg maximum

**System Height:**.....1677 mm (65.91 inches)  
For Floorstand

**Warranty Term:**.....2 Years Parts and Labor

#### Optional Accessories:

**Co-Observation:**.....Beam splitter with monocular  
observation available